

Features of *Vantage PT Software Version X3.1* by *PackTest Machines Inc.*

Sr No.	Feature	Description	Vantage PT Version X3.1
1	User Login	Unlimited Users with individually configurable privileges	Yes
		Incorrect Logins are logged. If user enters 3 wrong passwords, user is SUSPENDED and cannot login again even with correct password. Event is logged as Critical. Only an admin (or user with appropriate privileges can Un-Suspend the suspended user.	Yes
		Auto - Logout after 15 mins of inactivity. (Timeout period can be defined)	Yes
		Auto - Shutdown (software) on inactivity. (Auto-Logout will also trigger Auto-Shutdown)	Yes
2	User Rights / Privileges	Any user (with Modify user privileges) can create unlimited users with the any of following privileges.	
		1) Perform Test	Yes
		2) Search Test	Yes
		3) Calibrate Machine	Yes
		4) View Audit Trail	Yes
		5) Add / Edit / Delete users (User Management privilege)	Yes
		6) Configure Software	Yes
		7) Review and Approve Test / Calibration reports	Yes
3	User Levels	Unlimited Users with Unlimited variable privileges	Yes
		User List Print (summary or detailed)	Yes



4	Software Configuration Options	1) Change Test Folder	No
		2) Change Backup Folder	Yes
		3) Set Auto Backup On/Off	Yes
		4) Set Auto Backup Frequency	Yes
		5) Set Password Expiry On/Off	Yes
		6) Set Password Expiry Time (Days)	Yes
5	User Password Expiry / Change	Upon user first login, user has to change password immediately.	Yes
		User can change his/her own password whenever needed / expired.	Yes
		User has to change password after expiry date (Password expiry time can be defined by admin).	Yes
6	Password Complexity	Password should be 8-20 characters long and has to contain alphabets, numbers and &%#\$@.	Yes
		Complexity : Password should contain atleast 1 uppercase, 1 lower case a 1 special character	Yes
		New password cannot be same as previous 3 passwords.	Yes
7	User Deletion	User can be deleted by anyone who has User Management privilege.	Yes
		User cannot delete himself.	Yes
		No one can delete the SUPER USER.	Yes

8	User Suspension	User can be suspended/un-suspended by anyone who has Modify User privilege.	Yes
		Suspended user cannot login to software. If he tries to login, a Critical Audit Trail entry will be passed.	Yes
9	Data Storage	Data is stored in an internal database. PDF reports are generated as and when demanded by a user.	Yes
10	User Signature	User's details captured from the login is posted on the report "analysis done by" section.	Yes
11	Reviewer Signature	Reviewer can login with his credentials, review and approve Test and Calibration Data. Once approved, the report will carry his E-Sign.	Yes
12	Backup Data	If Autobackup is ON, software will automatically backup all configuration files and test data as per frequency chosen by IT Admin.	Yes
13	Audit Trail Data	Software captures important events and logs them in the audit trail log.	Yes
		A PDF report can be generated as and when required.	Yes
		Critical Events will be flagged separately.	Yes
14	Aborted Tests / Calibrations	Audit Trail Report can be pulled for a date range or user or event.	Yes
		Aborted Tests and Calibrations are saved and events logged in Audit Trail	Yes
15	Ignore a Test Value	No test result can be avoided from being reported. But user can mark a particular result as abnormal.	Yes
		Abnormal results are not considered in statistical calculations but will still show in report.	Yes
		User needs to justify the Abnormal values (if any) in the remarks section.	Yes
		A test report can have 10 normal samples + upto 5 abnormal samples. Only normal tests will be used for statistical calculations.	Yes
16	Test / Calibration Report Format	Report Format is fixed in software. Not changeable by user.	Yes

17	Report Header	Image Header with Logo or Text Header with 3 lines of Text.	Yes
18	Other Features	Incorrect Shut Downs are logged in Audit Trail.	Yes
		Equipment ID is printed on all Test / Calibration / Audit Trail Reports	Yes
		Automatic Pre-Calibration Diagnostics of communication.	Yes
		Option of Pre-Analysis Diagnostics (10-15 seconds) to prevent communication errors.	Yes
		Enhanced data handling and communication with equipment.	Yes